

Technical Specification Datasheet

Certificate of Conformity:

- The cell contained in this box meet the published spectral data for the raw material and the transmission characteristics specified in our published transmission chart.
- The cells have measured at 200nm, for T% > 80% and matched within 1.5%.
- Certified that the whole supplies detailed hereon have inspected tested and unless otherwise stated above, conform in all respects to the needs of the order description and drawings relative thereto. (Below. : Transmission certificate match pair)

TRANSMISSION CERTIFICATE			
CELL	% TRANSMISSION		
NO.	200nm	220nm	240nm
1.	81.8	82.5	84.1
2.	81.6	82.5	84.0

Techinstro assures the above specification and dimensions of the product as per datasheet. We ensure the quality of the product with moisture proof packing. This sheet is also available on the website, where a researcher or purchaser can also download it merely by clicking on the pdf link. The terms and condition for the product purchase mentioned in the quotation. For safe product storage, tips and procedure provided with the material. For more queries, contact our technical support.